Certified Ergonomics – For really relaxed work!



Transferpette® electronic

F I R S T C L A S S · B R A N D

Ergonomics – approved and certified.

The BRAND Transferpette® electronic pipettes combine the innovation of BRAND mechanical pipettes with electronic control and motorized operation to optimize ergonomics, accuracy and ease of use. All confirmed by the Technical Control Board of Rhineland with the Ergonomics Certificate.

Ergonomic

- Functional, ergonomic housing design
- Individually adjustable finger rest

Easy operation Intuitive menu structure

Innovative

Significantly reduced tip attachment and ejection forces

Convenient programs

- Pipetting
- Reverse pipetting
- Mixing
- GEL-Electrophoresis
- Dispensing
- Easy Calibration Technique

Ready for use

- 4000 pipetting cycles with one battery charge
- Battery refresh function

A	ERGONOMIE GEPRUFT
TÜV Rheinland Product Safety	ERGONOMICS APPROVED









The first worldwide!

The Transferpette[®] electronic was the first microliter pipette recognized with the ergonomics approved certificate from the Technical Control Board Rhineland/ Berlin-Brandenburg! Independent and neutral user tests confirm the ergonomics and the operating ease of the product and system! A user acceptance rating of **1.54** for the single channel Transferpette[®] electronic or **1.55** for the multichannel Transferpette[®] electronic is an outstanding result. You can obtain information about the single and multichannel Transferpette[®] electronic pipettes at www.tuv.com; ID No. 0011105500 and 5211207400.



Certified Ergonomics

Adapts to your hand

The optimal position of the thumb relative to the functional elements of the pipette is the starting point for a relaxed grip. The relaxed hand is an essential part of avoiding Repetitive Strain Injuries – RSI – from serial pipetting operations. The adjustable finger rest lets the Transferpette® electronic adapt to your hand for greatest comfort, whether you are leftor right-handed.

Easy handling – simple to operate

The menus are simple and obvious...there's no long learning curve or complicated programming. You just select the pipetting technique you've always used with mechanical pipettes, but benefit from a power-assist! The technical instructions are especially user friendly and always guide you through all functions with explanatory illustrations and straight-forward explanations.

Transferpette® electronic



Transferpette®-8 electronic



Lightweight

Innovative engineering and high quality materials preserve the light weight of the popular mechanical Transferpette[®] pipettes, while adding the control advantages of electronics.



Adjustable finger rest

Intensive and repeated work with poor ergonomic design can lead to muscular problems from the recurring stress. This is known as repetitive strain injury (RSI). In the laboratory, such injuries include tenosynovitis and carpal tunnel syndrome. To allow you to perform pipetting in a relaxed manner, the Transferpette[®] electronic has a continuously adjustable finger rest. This allows every user whether right- or left-handed, with large hands or small – to position the pipette in the hand so that the functional elements are within easy reach. The grip remains natural and comfortable for near effortless operation.

A Closer look..._


Models

Ergonomics - approved and certified.

Transferpette[®] electronic single channel pipette is available in 5 different models: 0.5-10 µl, 2-20 µl, 20-200 µl, 100-1000 µl and 0.5-5 ml.

The Transferpette®-8/-12 electronic multichannel pipette is available in 5 different models: 0.5-10 µl, 1-20 µl, 5-100 µl, 10-200 µl and 15-300 µl.

Designed for comfort in your relaxed hand, the Transferpette[®] electronic adds simple operating software and well-illustrated documentation to produce a uniquely user-friendly experience.





A Closer Look...

Large, clear display

Intuitive operation of all functions using 4 keys

Big pipetting button

Ergonomically arranged ejection button with color code for pipette volume and tip size

Transferpette®.8

No problem even with eight or twelve channels

The operation of 8- or 12-channels simultaneously multiplies the risk of RSI if the forces are not properly managed. The Transferpette® electronic manages them with a simple tap of the pipetting button; pipetting forces are uniform and negligible. A combination of innovations - including resilient V-rings on the tip cones, a stepped tip ejector, the manifold stabilizer and adjustable hand position all work together to creatively reduce ejection forces when working with a 8- or 12-channel Transferpette® electronic.



Individual shafts with seals can be easily unscrewed with only a simple gripping tool, which is supplied. Tip cones and seals can now be easily cleaned or replaced. This patented procedure eliminates the expense and long downtime associated with sending pipettes out for service, ensuring long service life and low operating costs. Recalibration is not required!

Complete manifold can be autoclaved at 121 °C and rotates freely 360° in either direction.

Stepped tip ejector reduces ejection forces significantly.

V-rings made of FKM allow effortless mounting and easy ejection of the tips, and readily adapt to tips from various manufacturers.

Tip cone with soft-seal components!



Two optional accessories are available for the single-channel Transferpette[®] models with volume ranges of 20-200 μ l and 100-1000 μ l. One is a patented shaft tip that is co-molded of hard and resilient components. Also available are three ejector clips These ensure that tips are not mounted with more force than needed to seal, so that tip ejection will take only a light touch.

The benefits of these accessories are especially noticeable in prolonged pipetting runs, where the cumulative strain of repetitive pipetting is the greatest risk.



Visual inspection! A resilient sealing ring not only adapts the pipette to small tip variations, but also provides a very distinct visual check for the correct seating and sealing of the tip.



Ejector Clips Not available on single channel models with volume ranges of 0.5 - 10 µl, 2 - 20 µl or 0.5 - 5 ml

Minimal Operating Force!

Freely rotatable over 360° in both directions! Complete manifold can be autoclaved!

The pipette manifold is freely rotatable with respect to the handle over 360° in both directions, so that you can always work at a relaxed and convenient angle to the titer plate.

Unscrew the manifold from the grip with a few quick turns, and autoclave the entire manifold at 121 °C without any further disassembly.



Stepped surface

FKM seal ring

The shaft and resilient, V-shaped FKM seal rings are designed so that only minimal attachment force is needed for solid and parallel tip seating. The stepped design allows the ejection force to be distributed to tips sequentially, but within fractions of a second, reducing ejection forces by 75% or more compared with unstepped ejector designs.

No Downtime!

The Transferpette[®] electronic draws its energy from an NiMH battery.

One battery charge allows you to perform over 4000 pipetting cycles.

During the charging process, you can continue to pipette, taking advantage of the handy one-meter charging cable.



Programs

Besides pipetting comfort, an important design goal for the Transferpette® electronic was rapid and intuitive work, and simple, fast program changes. The objective was to gain the advantages of electronic operation with no trade-offs. To simplify operation, complicated progams for rarely used functions were omitted. Thus, the Transferpette® electronic includes all of the functions used routinely in mechanical pipettes, while adding the ergonomic and performance advantages of power-assisted operation and digital accuracy. The functions of the Transferpette® electronic are simple to set and convenient to use.



Pipetting (PIP Mode)

The 'standard' program. The set volume is aspirated by the pipette, and then discharged.

Mixing of Samples (PIPmix Mode)

Program for mixing of liquids. The sample is repeatedly aspirated and discharged, and the number of mixing cycles is displayed.

Reverse Pipetting (revPIP Mode)

Program specially designed for the pipetting of liquids with a high viscosity, high vapor pressure or foamy media.

Pipetting with Electrophoresis (GEL Mode)*

Program for the loading of electrophoresis gels^{**}. The required sample volume is aspirated at the desired, adjustable speed, and is then discharged very slowly. The exact volume of liquid discharged is shown in the display as it is discharged.

Dispensing (DISP Mode)

A program for the dispensing of liquids in a series of equal aliquots. A volume that has been aspirated is dispensed in steps.

Easy Calibration Technique (CAL Mode)

Every laboratory that is certified according to ISO 9001 or accredited according to ISO/IEC 17025 or works according to the GLP/GMP guidelines must regularly calibrate its testing equipment and adjust as necessary. Program for making quick adjustments to the instrument, without tools. By changing the factory setting, <CAL> appears automatically in the display.

Battery-Refresh (batt Mode)

Regeneration function for increased performance and extending the service life of the batteries. The world's first microliter pipette with this function.


- The GEL mode is not included in the 1000 μl and 5000 μl single channel pipettes because these volumes are rarely used in electrophoresis.
- ** Patented













Tip Top Performance...



Pipette Tips and Filter Tips

Pipette tips, filter tips as well as primary packaging from BRAND are manufactured in a cleanroom under the most modern production conditions, and are automatically rack packed and packaged to ensure that the tips are of consistently high quality. For the manufacture, BRAND exclusively uses raw materials that are free from the additives di(2-hydroxyethyl)methyldodecylammonium (DiHEMDA) and 9-octadecenamide (oleamide).

- Standard- and Ultra Low Retention pipette tips
- High-purity polypropylene, free from DiHEMDA and oleamide
- Manufactured without lubricants
- Cadmium-free pigments
- Graduation for a quick volume check
- All racked tips and filter tips up to 1000 µl are free of DNA (<40 fg), RNase (<8.6 fg), endotoxins (<1 pg) and ATP (<1 fg)</p>
- Autoclavable at 121 °C (2 bar), acc. DIN EN 285
- Environmentally friendly packaging systems
- CE-marked according to IVD-Directive 98/79 EC
- Production is performed in Class 8 cleanroom, validated according to ISO 14 644-1.





Tip Top Performance...

TipStack™

TipRacks

TipBox - TipRack - TipStack[™] System

The tip system from BRAND has many distinctive features. The TipBox is equipped with transparent window for better all-round visual inspection. The side panel labeling on the colored tip-trays are always clearly visible.

The lid of the TipBox can be used as hinged and push-on lid. Simply rotate 180°.

Usable with a partially loaded multichannel pipette. The rimless tip-tray border enables problem-free loading of individual pipette tips with multichannel instruments. Tightly sealing spacers in the TipStack[™] prevent the tips from getting stuck together, and ensure them to be free from DNA, RNases, endotoxins and ATP. The sterile TipStacks and TipRacks are supplied with a transfer aid for contaminationfree use in a previously sterilized TipBox.

Cont

Transferpette® Transferpette® electronic electronic single channel multichannel Nominal volume Nominal volume Pipette tips Filter tips Volume range 크 Ξ Ξ 1000 µ п П 5000 Ξ 크 ᆿ ᆿ 200 200 100 300 9 20 9 20 0,1 - 20 µl 1 0,5 - 20 µl 1 - 50 µl 2 - 200 µl¹⁾ 5 - 300 µl 50 - 1000 µl¹⁾ 0,5 - 5 ml 1 - 10 ml 0,1 - 1 µl 0,5 - 10 µl 1 - 20 µl 2 - 20 ul 5 - 100 µl

Transfer aid

🖌 = Tip volume less than pipette's nominal volume

¹⁾ Bulk tips are yellow or blue colored, racked tips are colorless

in a yellow-colored or blue-colored tip tray

5 - 200 μl 50 - 1000 μl

10.0

²⁾ Calibration and adjustment may be necessary

Ordering date

Transferpette® electronic

Items supplied: each Transferpette[®] electronic pipette has DE-M marking and is supplied with performance certificate, battery, AC adapter, silicone grease.

Capacity, μl (color-coded)	Subdiv. µl	A* ≤ : %	± μl	CV* ≤ %	μl	With AC adapter for	Cat. No.
0.5 - 10	0.01	1.0	0.1	0.4	0.04	Europe (continental) UK/Ireland USA/Japan Australia without AC adapter	7052 99 7053 09 7053 19 7053 29 7053 39
• 2 - 20	0.02	1.0	0.2	0.4	0.08	Europe (continental) UK/Ireland USA/Japan Australia without AC adapter	7053 00 7053 10 7053 20 7053 30 7053 40
<mark>-</mark> 20 - 200	0.2	0.8	1.6	0.2	0.4	Europe (continental) UK/Ireland USA/Japan Australia without AC adapter	7053 03 7053 13 7053 23 7053 33 7053 43
• 100 - 1000	1.0	0.6	6	0.2	2	Europe (continental) UK/Ireland USA/Japan Australia without AC adapter	7053 06 7053 16 7053 26 7053 36 7053 46
5 00 - 5000	5.0	0.6	30	0.2	10	Europe (continental) UK/Ireland USA/Japan Australia without AC adapter	7053 07 7053 17 7053 27 7053 37 7053 47

Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are well within the limits of DIN EN ISO 8655-2. DE-M marking. A = Accuracy, CV = Coefficient of variation

Note! When ordering instruments with DAkkS calibration certificates, the prefix 'DAkkS' must be added to the order number, e.g., DAkkS 7052 99.

BRAND also offers calibration service at the factory lab.

Ordering data

Transferpette®-8 electronic

Items supplied: each Transferpette[®]-8 electronic pipette has DE-M marking and is supplied with performance certificate, battery, AC adapter, device stand, TipBox, refill unit, reagent reservoir, mounting tool, silicone grease and 1 set of sealing rings made of FKM.

Capacity, µl (color-coded)	Subdivision µl	A* ≤ ± % µl	CV* ≤ % µl	With AC adapter for	Cat. No.
0.5 - 10	0.01	1.2 0.12	0.8 0.08	Europe (continental) UK/Ireland USA/Japan	7053 99 7054 09 7054 19
1 - 20	0.02	1.0 0.2	0.5 0.1	Australia Europe (continental) UK/Ireland USA/Japan Australia	7054 29 7054 00 7054 10 7054 20 7054 30
<mark>9</mark> 5 - 100	0.1	0.8 0.8	0.25 0.25	Europe (continental) UK/Ireland USA/Japan Australia	7054 03 7054 13 7054 23 7054 33
<mark>-</mark> 10 - 200	0.2	0.8 1.6	0.25 0.5	Europe (continental) UK/Ireland USA/Japan Australia	7054 04 7054 14 7054 24 7054 34
• 15 - 300	0.5	0.6 1.8	0.25 0.75	Europe (continental) UK/Ireland USA/Japan Australia	7054 06 7054 16 7054 26 7054 36



Transferpette®-12 electronic

Items supplied: each Transferpette[®]-12 electronic pipette has DE-M marking and is supplied with performance certificate, battery, AC adapter, device stand, TipBox, refill unit, reagent reservoir, mounting tool, silicone grease and 1 set of sealing rings made of FKM.

Capacity, µl (color-coded)	Subdivision µl	A* ≤ ± % µl	CV* ≤ % µl	With AC adapter for	Cat. No.
0.5 - 10	0.01	1.2 0.12	0.8 0.08	Europe (continental) UK/Ireland USA/Japan Australia	7054 49 7054 59 7054 69 7054 79
1 - 20	0.02	1.0 0.2	0.5 0.1	Europe (continental) UK/Ireland USA/Japan Australia	7054 50 7054 60 7054 70 7054 80
<mark>-</mark> 5 - 100	0.1	0.8 0.8	0.25 0.25	Europe (continental) UK/Ireland USA/Japan Australia	7054 53 7054 63 7054 73 7054 83
0 10 - 200	0.2	0.8 1.6	0.25 0.5	Europe (continental) UK/Ireland USA/Japan Australia	7054 54 7054 64 7054 74 7054 84
15 - 300	0.5	0.6 1.8	0.25 0.75	Europe (continental) UK/Ireland USA/Japan Australia	7054 56 7054 66 7054 76 7054 86



* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are well within the limits of DIN EN ISO 8655-2. DE-M marking. A = Accuracy, CV = Coefficient of variation

Accessories

3-device charging stand for Transferpette[®] electronic (up to 1000 μl)

You can use the 3-device charging stand to charge up to three single channel pipettes at the same time. Pack of 1.

for Transferpette [®] electronic with AC adapter for	Cat. No.	
Europe (continental) (230V/50Hz)	7053 90	
UK/Ireland (230V/50Hz)	7053 91	

UK/Ireland (230V/50Hz)	7053 91
USA/Japan (110V/50-60	Hz) 7053 92
Australia (240V/50Hz)	7053 93




Individual stand for Transferpette[®] electronic

Also suitable for corresponding models of Transferpette[®] and Transferpette[®] *S* pipettes. Pack of 1.

for Transferpette [®] electronic	Cat. No.
up to 1000 µl	7053 85
500-5000 µl	7053 86







Cat. No.	7034 40

For all Transferpette® multi-

Individual stand

channel pipettes. Pack of 1

Reagent reservoir, PP

Non-sterile, with lid. Pack of 10.


Cat. No. 7034 59

Sterile, without lid. Packed individually. Pack of 100.

Cat. No.	7034 11
Sterile, without I	id.
5 ner hag	

7034 09

5 per bag. Pack of 200. Cat. No.



Accessories for Transferpette[®] electronic 20-200 µl and 100-1000 µl

Description	Cat No.
 Ejector-adjustment clips (set of 3) 20 - 200 μl Ejector-adjustment clips (set of 3) 100 - 1000 μl 	7055 60 7055 62
 Pipette shaft (complete) 20 - 200 μl Pipette shaft (complete) 100 - 1000 μl 	7055 43 7055 45

 $\mathsf{BRAND}^{\otimes},\,\mathsf{BIO}\text{-}\mathsf{CERT}^{\otimes}$ and Transferpette^ are trademarks of BRAND GMBH + CO KG, Germany. Other reproduced brands are the property of the respective owner.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

Subject to technical modification without notice. Errors excepted.

